

Nalidixic acid

CMS213


Product Identifier

CAS No.	:	389-08-2
EC No.	:	206-864-7
Molecular Formula	:	C ₁₂ H ₁₂ N ₂ O ₃
Molecular Weight	:	232.24
HS Code	:	2933 99 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow to yellow-green crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of dichloromethane
FTIR	:	Matches with the standard pattern
Melting range	:	224 - 233°C
Assay (NaOH Titration/HPLC)	:	>= 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H351
Precautionary Statement(s)	:	P270- P280- P301 + P312- P308 + P313- P330- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QN2885000
WGK	:	1