

## Menadione sodium bisulphite, water soluble

**CMS3031**

### Product Identifier

CAS No.	:	130-37-0
EC No.	:	204-987-0
EC Index No	:	607-618-00-5
Molecular Formula	:	C <sub>11</sub> H <sub>9</sub> NaO <sub>5</sub> S
Molecular Weight	:	276.24
Synonym	:	Vitamin K <sub>3</sub> sodium bisulfite, water soluble
HS Code	:	2914 79 90
Storage	:	2 - 8°C
Shelf life	:	3 years

### Technical Specification

Appearance	:	White to yellow hygroscopic crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
Water (K.F.)	:	<= 15%
Assay (HPLC)	:	>= 95.00%

### Safety Information

Hazard Pictogram(s) :



Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H410
Precautionary Statement(s)	:	P273- P280- P337 + P313- P391- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	Not available
WGK	:	2