

Cefaclor

CMS3310


Product Identifier

CAS No.	:	53994-73-3
EC No.	:	258-909-5
Molecular Formula	:	C ₁₅ H ₁₄ ClN ₃ O ₄ S
Molecular Weight	:	367.81
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to light yellow powder
Solubility	:	33.3 mg soluble in 1 mL of 1N hydrochloric acid
FTIR	:	Matches with the standard pattern
Assay (HPLC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317
Precautionary Statement(s)	:	P261- P280- P321- P333 + P313- P362 + P364- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XI0363000
WGK	:	2