

## Oxacillin sodium salt monohydrate

**CMS5372**


### Product Identifier

CAS No.	:	7240-38-2
EC No.	:	214-636-3
Molecular Formula	:	C <sub>19</sub> H <sub>18</sub> N <sub>3</sub> O <sub>5</sub> SNa.H <sub>2</sub> O
Molecular Weight	:	441.43
HS Code	:	2941 10 20
Storage	:	2 - 8°C
Shelf life	:	3 years

### Technical Specification

Appearance	:	White to almost white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (3% in water at 25°C)	:	4.50 - 7.50
Potency	:	815 - 950 µg/mg oxacillin

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H317- H319- H334- H335
Precautionary Statement(s)	:	P261- P280- P305+P351+P338- P342 + P311
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3