

Novobiocin sodium salt

CMS643


Product Identifier

CAS No.	:	1476-53-5
EC No.	:	216-023-6
Molecular Formula	:	$C_{31}H_{35}N_2NaO_{11}$
Molecular Weight	:	634.61
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow hygroscopic crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Potency	:	$\geq 850 \mu\text{g}/\text{mg}$

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H317- H319
Precautionary Statement(s)	:	P261- P280- P333 + P313- P337 + P313- P362 + P364- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RD5425000
WGK	:	1