

CMS7464

Rutin trihydrate, For Biochemistry

Product Number

Packing

CMS7464 : 100G
CMS7464 : 25G

Product Information

Product Code : CMS7464
Product Name : Rutin trihydrate, For Biochemistry
Synonym : Vitamin P trihydrate
Molecular Formula : $C_{27}H_{30}O_{16} \cdot 3H_2O$
Molecular Weight : 664.56
CAS No. : 250249-75-3
EC No. : 205-814-1
HS Code : 2938 10 00
Shelf Life : 3 years

Technical Specification

Appearance : Off-white to yellow to green crystals or powder
Solubility : Insoluble in water
Assay (TBAH Titration) : 90.00 - 101.00%


GHS Safety Information

Hazard Statement(s) : H302
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 22
WGK : 3
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods