

Puromycin dihydrochloride, From Streptomyces alboniger

CMS8861


Product Identifier

CAS No.	:	58-58-2
EC No.	:	200-387-8
Molecular Formula	:	C ₂₂ H ₂₉ N ₇ O ₅ .2HCl
Molecular Weight	:	544.44
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to beige crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (HPLC)	:	min. 98.50%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302
Precautionary Statement(s)	:	P301+P312+P330
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AU7355000