

Sulfathiozole

CMS8929


Product Identifier

CAS No.	:	72-14-0
EC No.	:	200-771-5
Molecular Formula	:	C ₉ H ₉ N ₃ O ₂ S ₂
Molecular Weight	:	255.32
HS Code	:	2935 00 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to cream crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of 1M NaOH
FTIR	:	Matches with the standard pattern.
Melting range	:	199 - 202°C
Assay (NT/NaNO ₂ Titration)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WP2360000
WGK	:	2