
Teicoplanin

CMS8940

Product Identifier

CAS No. : 61036-62-2
HS Code : 3822 00 90
Storage : 2 - 8°C
Shelf life : 3 years

Technical Specification

Appearance : White to yellow powder or solid
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 6.30 - 7.70
Potency : ≥ 900 $\mu\text{g}/\text{mg}$

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : WY1559000
WGK : 1