

Timentin

CMS9138

Ticarcillin disodium salt : Potassium clavulanate mixture (15:1)

Product Identifier

| | | |
|------------|---|------------|
| CAS No. | : | 86482-18-0 |
| HS Code | : | 2941 90 90 |
| Storage | : | 2 - 8°C |
| Shelf life | : | 3 years |

Technical Specification

| | | |
|-------------------------------|---|----------------------------------|
| Appearance | : | White to yellow powder |
| Solubility | : | 33.3 mg soluble in 1 mL of water |
| pH (0.005 % in water at 25°C) | : | 5.50 - 7.50 |

Safety Information

Hazard Pictogram(s)



| | | |
|----------------------------|---|---|
| Signal Word | : | Danger |
| Hazard Statement(s) | : | H315- H319- H334- H335 |
| Precautionary Statement(s) | : | P261- P280- P284- P304 + P340- P312- P501 |
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 1 |