

## Clindamycin hydrochloride

**CMS9386**


### Product Identifier

CAS No.	:	21462-39-5
EC No.	:	244-398-6
Molecular Formula	:	C <sub>18</sub> H <sub>33</sub> ClN <sub>2</sub> O <sub>5</sub> S.HCl
Molecular Weight	:	461.44
HS Code	:	2941 90 90
Storage	:	Below 30°C
Shelf life	:	3 years

### Technical Specification

Appearance	:	White to almost white crystals or powder
Solubility	:	100 mg soluble 1 mL of water
pH (10% in water at 25°C)	:	3.00 - 5.50
Assay (HPLC)	:	>= 96.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H317- H319- H362
Precautionary Statement(s)	:	P260- P263- P280- P308 + P313- P333 + P313- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1