

DEV Lactose Peptone Broth, Granulated

GM1356

DEV Lactose Peptone Broth, granulated is recommended for enrichment and titre determination of coliform bacteria.

Composition**

Ingredients	Gms / Litre
Meat peptone	10.000
Lactose	10.000
Sodium chloride	5.000
Bromo cresol purple	0.010
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 25.01 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Dispense in tubes or flasks as desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle And Interpretation

DEV Lactose Peptone Broth is used for enrichment and titre determination of coliform bacteria.

Meat peptone provides carbon and other essential nutrients for growth. Sodium chloride is for the osmotic balance Lactose is the carbohydrate source. Bromocresol purple is a pH indicator which has a yellow colour below pH 5.3 and purple colour above pH 6.7. Due to the fermentation of lactose and the subsequent acid production the colour changes to yellow (1).

Quality Control

Appearance

Cream to yellow coloured granular medium

Colour and Clarity of prepared medium

Purple coloured clear solution

Reaction

Reaction of 2.5 % w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

Cultural characteristics observed after an incubation at 37°C for 18-24 hours.

Organism

Growth

Cultural Response

Enterobacter aerogenes ATCC 13048 good

Escherichia coli ATCC 25922 good

Klebsiella pneumoniae ATCC 13883 good

Storage and Shelf Life

Store below 30°C in tightly closed container and use freshly prepared medium. Use before expiry date on the label.

Reference

1.MacFaddin J.F., 1985, Media for Isolation - Cultivation - Identification - Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

Revision : 00 / 2014



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.