

Gelatin DEV Agar, Granulated

GM1609

Gelatin DEV Agar, granulated is used for determining the total microbial count and detecting gelatin-liquefying microorganisms in water.

Composition**

Ingredients	Gms / Litre
Peptone from meat	10.000
Meat extract	10.000
Sodium chloride	5.000
Gelatin	10.000
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 50 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 121°C for 15 minutes. DO NOT OVERHEAT. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Gelatin, a protein derivative of animal collagen is incorporated into various media to determine an organisms ability to produce proteolytic type enzyme (proteinase) detected by digestion or liquefaction of the gelatin (1). Gelatinase is a pepsin, which hydrolyses gelatin to polypeptides, peptides and amino acids. Gelatin DEV Agar is used to determine the total microbial count and for detecting gelatin liquefying microorganisms in water as per the German Drinking Water Regulations, 1990 (2).

The medium consists of nutrients like peptone from meat, meat extract and gelatin, which provide nitrogen compounds and also carbon compounds for the growth of organisms. Gelatin acts as solidifying agent and also acts as a substrate for the organisms producing gelatinase enzyme. Gelatin breakdown can be visualized by flooding the plates with a saturated solution of ammonium sulphate. Clear zones are observed around gelatin-liquefying colonies.

Quality Control

Appearance

Cream to yellow coloured granular medium

Gelling

Firm, comparable with 1.5% Agar gel and 1.0% Gelatin gel

Colour and Clarity of prepared medium

Yellowish brown coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 5.0% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH

7.10-7.50

Cultural Response

Cultural characteristics observed after an incubation at 18-22°C for 40-48 hours.(Gelatin liquefaction is observed by flooding the plate with saturated solution of ammonium sulphate)

Cultural Response

Organism	Inoculum (CFU)	Growth	Recovery	Gelatin liquefaction
Cultural Response <i>Escherichia coli</i> ATCC 25922	50-100	good-luxuriant	≥70%	negative, no clear zone

<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	good-luxuriant	$\geq 70\%$	around the colony positive, clear zone around the colony
<i>Bacillus cereus</i> ATCC 10876	50-100	good-luxuriant	$\geq 70\%$	positive, clear zone around the colony
<i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant	$\geq 70\%$	positive, clear zone around the colony
<i>Proteus vulgaris</i> ATCC 13315	50-100	good-luxuriant	$\geq 70\%$	negative, no clear zone around the colony
<i>Enterococcus faecalis</i> ATCC 29212	50-100	good-luxuriant	$\geq 70\%$	positive, clear zone around the colony
<i>Aeromonas hydrophila</i> ATCC 7966	50-100	good-luxuriant	$\geq 70\%$	positive, clear zone around the colony

Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

Reference

1. MacFaddin J. F., 2000, Biochemical tests for Identification of Medical Bacteria, 3rd Edi., Lippincott, Williams and Wilkins, Baltimore.
2. Deutsche Einheitsverfahren zur Wasser-, Abwasser- und Schlammuntersuchung. - VCH Verlagsgesellschaft, D-6940 Weinheim.
3. Verordnung über Trinkwasser und über Wasser für Lebensmittelbetriebe vom 12, Dezember 1990, - Bundesgesetzbl.: Teil I; 2613-2669 (1990).

Revision : 00 / 2014



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.