

Alicyclobacillus Agar, Granulated

GM1650

Alicyclobacillus Agar, granulated is used for the isolation of *Alicyclobacillus* species in fruit juices.

Composition**

Ingredients	Gms / Litre
Yeast extract	2.000
Dextrose	5.000
Calcium chloride	0.25066
Magnesium sulphate	0.500
Ammonium sulphate	0.200
Potassium dihydrogen phosphate	3.000
Zinc sulphate	0.00018
Copper sulphate	0.00016
Manganese sulphate	0.00015
Cobalt chloride	0.00018
Boric acid	0.00010
Sodium molybdate	0.00030
Agar	18.000
Final pH (at 25°C)	4.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 28.95 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Note: Adjust the pH of the medium to 4.0 ± 0.2 (after sterilization) using 1N H₂SO₄ or 1N NaOH .

Principle And Interpretation

Alicyclobacillus species are gram positive aerobic thermophilic, and spore forming acidophilic bacteria. *Alicyclobacillus* are sometimes called Acidophilic Thermophilic Bacteria (ATB). These spore forming organisms are able to survive the relatively mild pasteurization temperatures used for fruit juices and drinks and some are able to grow out and cause spoilage of the beverage. Even very low numbers of *Alicyclobacillus* are able to cause spoilage and produce objectionable flavours and odours specially affecting the quality of fruit juice (1,2) and in the beverages, damaging the brand. These bacteria are able to grow at pH values as low as 2.5 and also at elevated temperatures as high as 60°C.

Alicyclobacillus Agar has a pH of 4.0 ± 0.2 which inhibits most of the microbial flora. Rest of the microbial flora is inhibited at 60°C, which is the optimum growth temperature for *Alicyclobacillus* species. Alicyclobacillus Agar is a slight modification of the medium recommended for the cultivation of *Alicyclobacillus acidoterrestris* (3).

Quality Control

Appearance

Cream to yellow coloured granular medium

Gelling

Firm, comparable with 1.8% Agar gel

Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 2.89% w/v aqueous solution at 25°C. pH : 4.0±0.2

pH

3.80-4.20

Cultural Response

Cultural characteristics observed after an incubation at 60°C for 48-72 hours.

Cultural Response

Organism	Growth
Cultural Response	
<i>Alicyclobacillus</i>	good to
<i>acidocaldarius</i> 27009	luxuriant
<i>Escherichia coli</i> ATCC 25922	inhibited
<i>Staphylococcus aureus</i> ATCC 25923	inhibited
<i>Saccharomyces cerevisiae</i> ATCC 19615	inhibited
<i>Candida albicans</i> ATCC 10231	inhibited

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label.

Reference

1. Ceny G., Hennlich W. and K Rocallia-Furchtsaftwerb duchr Baciilen. Isobioerung and Charakteriseeuing des Verdebserregens-Z hebers Utres Forsch 179: 224-227, 1984.
2. Baungart and Merve S., The Impact of Alicyclobacillus acidoterstris on the Quality of Juices and Soft Drinks Fruit processing 7: 251-254 (2000).
3. Atlas R. M., 2010, Handbook of Microbiological Media, 4th Ed, CRC Press.

Revision : 00 / 2014



Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.