

## Giolitti Cantoni Broth Base, Granulated

**GM584**

Giolitti Cantoni Broth Base, granulated is used for selective enrichment of *Staphylococcus aureus* from foods.

### Composition\*\*

<b>Ingredients</b>	<b>Gms / Litre</b>
Casein enzymic hydrolysate	10.000
Beef extract	5.000
Yeast extract	5.000
Mannitol	20.000
Sodium chloride	5.000
Lithium chloride	5.000
Glycine	1.200
Sodium pyruvate	3.000
Final pH ( at 25°C)	6.9±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 54.2 grams in 1000 ml distilled water. Warm gently to dissolve the medium completely. Dispense 19 ml amounts in 20mm x 200mm test tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool rapidly to room temperature and aseptically add 0.3 ml of 3.5% Potassium Tellurite Solution (FD047) to each tube. Add 0.03 ml for testing meat and meat products. Mix well before use.

*Caution : Lithium chloride is harmful. Avoid bodily contact and inhalation of vapours. On contact with skin, wash with plenty of water immediately .*

### Principle And Interpretation

Giolitti-Cantoni Broth Base is a fluid medium employed for the recovery of low number of Staphylococci from foodstuffs as described by Giolitti and Cantoni (6). Giolitti- Cantoni Broth was also recommended by Mossel et al for detecting *Staphylococcus aureus* in dried milk, baby food and other food products (1). This medium was recommended as an enrichment medium by the International Dairy Federation (IDF) and APHA for detecting *S.aureus* in dried milk and other foods stating that the organism should be absent in 1 gram of sample (2,3). ISO committee has also recommended this medium for examination of meat and meat products (4).

Giolitti-Cantoni Broth Base contains casein enzymic hydrolysate, yeast extract and beef extract as sources of carbon, nitrogen, vitamins and minerals and B-complex vitamins. Mannitol and sodium pyruvate in the basal medium act as growth stimulants for *S.aureus* . Lithium chloride inhibits gram-negative lactose fermenting bacilli. Potassium tellurite and glycine inhibit gram-positive bacilli (7). Addition of sterile paraffin wax to the inoculated medium inhibits Micrococci due to creation of anaerobic conditions. Potassium tellurite concentration must be reduced as per the weight of test sample (0.1 - 0.01 gram). The medium should be inoculated as soon as it has been cooled after sterilization, otherwise absorbed oxygen should be expelled by placing the tubes in free-flowing steam for 15-20 minutes.

Inoculate 1 gram of sample or 1 ml of a suitable dilution of a sample into 19 ml of Giolitti-Cantoni Broth tubes in duplicate. Overlay the medium with 5 ml molten sterile paraffin wax and incubate at 37°C for 24-48 hours and examine daily. Blackening of the medium (usually at the bottom) within 48 hours indicates the presence of *S.aureus*. The blackened medium, when streaked on Baird Parker Agar (GM043/M043), shows black colonies surrounded by clear zones (5).

### Quality Control

#### Appearance

Light yellow to brownish yellow coloured granular medium

#### Colour and Clarity of prepared medium

Medium amber coloured, clear solution without any precipitate

### Reaction

Reaction of 5.42% w/v aqueous solution at 25°C. pH : 6.9±0.2

### pH

6.70-7.10

### Cultural Response

Cultural characteristics observed with added 3.5% Potassium Tellurite Solution (FD047), after an incubation at 35-37°C for 24-48 hours.

### Cultural Response

Organism	Inoculum (CFU)	Growth	Tellurite reduction
<b>Cultural Response</b> <i>Staphylococcus aureus</i> ATCC 25923	50-100	good-luxuriant	positive, blackening at the bottom of the tubes or general blackening of the medium
<i>Escherichia coli</i> ATCC 25922	≥10 <sup>3</sup>	inhibited	
<i>Micrococcus luteus</i> ATCC 10240	≥10 <sup>3</sup>	inhibited	
<i>Staphylococcus aureus</i> ATCC 6538	50-100	good-luxuriant	positive, blackening at the bottom of the tubes or general blackening of the medium
<i>Bacillus cereus</i> ATCC 11778	≥10 <sup>3</sup>	inhibited	
<i>Pseudomonas aeruginosa</i> ATCC 27853	≥10 <sup>3</sup>	inhibited	
<i>Staphylococcus epidermidis</i> ATCC 12228	50-100	poor-fair	variable reaction

### Storage and Shelf Life

Store below 30°C in tightly closed container and prepared medium at 2-8°C. Use before expiry period on the label.

### Reference

1. Mossel D. A. A., Harrewijn G. A. and Elzebroek J. M., 1973, UNICEF.
2. International Dairy Federation, 1978, IDF Standard 60A:1978, International Dairy Federation, Brussels, Belgium.
3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
4. International Organization for Standardization (ISO), 1977, Draft ISO/DIS 5551, Part 2.
5. De Waart J., Mossel D. A. A., Ten Broeke R. and Van de Moosdijk A., 1968, J. Appl. Bacteriol., 31:276.
6. Giolitti C. and Cantoni C., 1966, J. Appl. Bacteriol., 29: 395.
7. Lambin S. and German A., 1961, Precis De Microbiologie, pg. 63, Paris Masson.

Revision :00 / 2014



### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.