
Sodium chloride, Hi-LR™

GRM031

Product Identifier

CAS No.	:	7647-14-5
EC No.	:	231-598-3
Molecular Formula	:	NaCl
Molecular Weight	:	58.44
HS Code	:	2501 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white crystals or powder or solid
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	5.00 - 9.00
Sulphate (SO ₄)	:	<= 0.02%
Iron (Fe)	:	<= 0.002%
Assay (AT)	:	98.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VZ4725000
WGK	:	1