

---

## Papain, Purified, From *Papaya Latex*

**GRM058**

### Product Identifier

---

CAS No.	:	9001-73-4
EC No.	:	232-627-2
EC Index No	:	647-007-00-0
HS Code	:	3507 90 71
Storage	:	2 - 8°C
Shelf life	:	3 years


### Technical Specification

---

Appearance	:	White to light brown amorphous powder
Solubility	:	10 mg soluble in 1 mL of water
Assay (Proteolytic activity)	:	The difference between Test & Blank titration is not less than 4.5 ml of 0.1 M NaOH

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H319- H334- H335
Precautionary Statement(s)	:	P261- P280-P304 + P340- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RU4950000
WGK	:	1