
Albumin egg, flakes

GRM080

Product Identifier

CAS No.	:	9006-59-1
EC No.	:	232-692-7
HS Code	:	3502 19 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow to dark brown flakes
Solubility	:	In water, it swells at first and is sparingly soluble
Loss on drying (at 105°C, 1 hr)	:	<= 15.00%
Sulphated ash	:	<= 8.00%
Total nitrogen	:	13.00 - 16.00%
Protein content (on dry basis)	:	>= 80.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H334
Precautionary Statement(s)	:	P261- P284- P304 + P340- P342 + P311- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AY9384000
WGK	:	1