
Potassium tellurite, Hi-LR[™] (Suitable for bacteriological work)

GRM090

Product Identifier

CAS No.	:	7790-58-1
EC No.	:	232-213-1
Molecular Formula	:	K ₂ TeO ₃
Molecular Weight	:	253.79
HS Code	:	2842 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to grayish white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water to an alkaline solution
Assay (Gravimetry)	:	min. 90.00%

Safety Information

Hazard Pictogram(s) :



Signal Word	:	Danger
Hazard Statement(s)	:	H301- H315- H319- H335
Precautionary Statement(s)	:	P280- P301+P310+P330- P304+P340+P312- P305+P351+P338- P337 + P313
UN No.	:	3284
Class	:	6.1
Packing Group	:	III
RTECS	:	Not available
WGK	:	3