

Urease

GRM091


Product Identifier

CAS No.	:	9002-13-5
EC No.	:	232-656-0
Synonym	:	Jackbean meal
HS Code	:	3507 90 99
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Cream to brown coloured powder or flakes or solid
Solubility	:	Insoluble in water
Activity	:	Each mg of jackbean meal (urease active meal) will hydrolyse 3 mg ureain 30 minutes at 37°C

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H319- H334- H335
Precautionary Statement(s)	:	P261- P280- P304 + P340- P342 + P311- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	YU1700000
WGK	:	1