
Xylenol orange tetrasodium salt, Practical grade

GRM1006

Product Identifier

CAS No.	:	3618-43-7
EC No.	:	222-805-8
Molecular Formula	:	$C_{31}H_{28}N_2O_{13}SNa_4$
Molecular Weight	:	760.58
HS Code	:	2934 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Dark red to dark brown crystals or powder
Solubility	:	10 ppm in 0.1 N NaOH yields clear solution
Absorption maximum	:	575 - 585 nm (10 ppm in 0.1 M NaOH)
Loss on drying (at 110°C, 1 hr)	:	<= 10.00%
Suitability for Zinc titration	:	Passes test

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1