

2,4-Dinitroaniline

GRM10168


Product Identifier

CAS No.	:	97-02-9
EC No.	:	202-553-5
Molecular Formula	:	C ₆ H ₅ N ₃ O ₄
Molecular Weight	:	183.12
HS Code	:	2921 42 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to gold to brown crystals or powder or chunks
Solubility	:	10 mg soluble in 1 mL of methanol
FTIR	:	Matches with the standard pattern
Melting range	:	176 - 182°C
Assay (GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300- H310- H330- H373- H411
Precautionary Statement(s)	:	P260- P264- P273- P280- P284- P301 + P310
UN No.	:	1596
Class	:	6.1
Packing Group	:	II
RTECS	:	BX9100000
WGK	:	2