

GRM10186 Ethyl cellosolve®**Product Number Packing**

GRM10186 : 500ML

Product Information

Product Code : GRM10186
Product Name : Ethyl cellosolve®
Synonym : Ethylene glycol monoethyl ether; 2-Ethoxyethanol
Molecular Formula : C₄H₁₀O₂
Molecular Weight : 90.12
CAS No. : 110-80-5
EC No. : 203-804-1
HS Code : 2909 44 00
Other Information : Key: (R) = Registered trademark of Union Carbide Corp.
Shelf Life : 4 years

Technical Specification

Appearance : Colourless liquid
Solubility : 33.3 mg soluble in 1 mL of water
FTIR (Liquid film) : Matches with the standard pattern
Refractive index (n_{20/D}) : 1.4060 - 1.4100
Density (at 25°C) : 0.925 - 0.935 g/mL
Assay (GC) : min. 99.00%

GHS Safety Information

Hazard Statement(s) : H226-H302-H312-H332-H360FD
Precautionary Statement(s) : P201-P280-P308+P313
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 60-61-10-20/21/22
S-Phrase(s) : 53-45
WGK : 2
RTECS : KK8050000
Flash Point(°F) : 107.6 °F
Flash Point(°C) : 42 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)





Product Information

Transport Information

UN No.	:	1171
Class	:	3
Packaging Group	:	3
Marine Pollutant	:	No
ADR/RID	:	1171 3/PG 3
IMDG	:	1171 3/PG 3
IATA	:	1171 3/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein