
Inulin, Hi-AR[™]

GRM103

Product Identifier

CAS No.	:	9005-80-5
EC No.	:	232-684-3
Molecular Formula	:	(C ₆ H ₁₀ O ₅) _n
HS Code	:	1108 20 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of hot water
pH (10% in hot water at 25°C)	:	4.00 - 7.00
Loss on drying (at 105°C, 2 hr)	:	<= 10.00%
Specific rotation	:	-37° to -15° (c = 1% in hot water, at 20°C)

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1