
Paraplast® For general tissue embedding

GRM10301

Product Identifier

CAS No.	:	8002-74-2
EC No.	:	232-315-6
HS Code	:	2712 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white pellets or solid.
Solubility	:	Insoluble in water.
Melting Point	:	~56°C

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RV0350000
WGK	:	nwg