

---

## Poly(propylene glycol) MW 2000

**GRM10321**

### Product Identifier

---

CAS No.	:	25322-69-4
EC No.	:	500-039-8
Molecular Formula	:	$H[OCH(CH_3)CH_2]_nOH$
Molecular Weight	:	~ 2000
HS Code	:	3907 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colorless viscous liquid
Solubility	:	1 mL miscible in 1 mL of methanol
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 0.10%
Acid value	:	<= 0.1 mg KOH/g
Hydroxyl value	:	53.00 - 58.00 mg KOH/g

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	TR5250000
WGK	:	1