
Sunset yellow

GRM10357


Product Identifier

CAS No.	:	2783-94-0
EC No.	:	220-491-7
C.I. No.	:	15985
Molecular Formula	:	$C_{16}H_{10}N_2Na_2O_7S_2$
Molecular Weight	:	452.37
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Orange-red to red to very dark red powder or granules or crystals
Solubility	:	10 ppm in water yields clear solution
Absorption maximum (10 ppm in water)	:	477 - 487 nm
Loss on drying (at 105°C, 2 hr)	:	<= 10.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QK2450000
WGK	:	2