

**GRM10402 Buffer solution (Borate), pH 10.00 ± 0.05****Product Number Packing**

GRM10402 : 1L

**Product Information**

Product Code : GRM10402  
Product Name : Buffer solution (Borate), pH10.00 ± 0.05  
Synonym : T.T.B  
HS Code : 3822 00 00  
Shelf Life : 3 years

**Technical Specification**

Appearance : Clear colourless liquid  
Solubility : 1 mL miscible in 1 mL of water  
pH (at 25°C) : 9.95 - 10.05

**GHS Safety Information**

Hazard Statement(s) : H315-H319-H360  
Precautionary Statement(s) : P201-P305 + P351 + P338-P308 + P313  
Signal Word : Danger  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/38  
S-Phrase(s) : 26  
WGK : 1  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods