

---

---

## Paraffin wax (noncaking), Packed in SS container

**GRM1042**

### Product Identifier

---

CAS No. : 8002-74-2  
EC No. : 232-315-6  
HS Code : 2712 20 00  
Storage : Below 30°C  
Shelf life : 4 years

### Technical Specification

---

Appearance : Colourless or white translucent mass  
Solubility : 33.3 mg soluble in 1 mL of chloroform  
Congealing Point : 60 - 62°C

### Safety Information

---

UN No. : Not dangerous goods  
Class : -  
Packing Group : -  
RTECS : RV0350000  
WGK : nwg