

**GRM10460 Cesium chloride, Hi-LR™**

**Product Number Packing**

---

GRM10460 : 100G  
GRM10460 : 500G

**Product Information**

---

Product Code : GRM10460  
Product Name : Cesium chloride, Hi-LR™  
Molecular Formula : CsCl  
Molecular Weight : 168.36  
CAS No. : 7647-17-8  
EC No. : 231-600-2  
HS Code : 2827 39 85  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : Colourless or white deliquescent crystals or powder or beads or chunks  
Solubility : 1000 mg soluble in 1 mL of water  
Assay (AT/EDTA Titration) : 99.00 - 102.00%

**Risk and Safety Information**

---

WGK : 2  
RTECS : FK9625000  
Storage Temperature(°C) : Store below 30°C

**Transport Information**

---

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods