

RM10574 Isoamyl acetate

Product Number Packing

RM10574 : 500ML

Product Information

Product Code : RM10574
 Product Name : Isoamyl acetate
 Synonym : Isopentyl acetate; Acetic acid 3-methylbutyl ester
 Molecular Formula : C₇ H₁₄ O₂
 Molecular Weight : 130.18
 CAS No. : 123-92-2
 EC No. : 204-662-3
 HS Code : 2915 39 00
 Shelf Life : 4 years

Technical Specification

Appearance : Colourless liquid
 Solubility : 33.3 mg soluble in 1 mL of ethanol
 FTIR (Liquid film) : Matches with the standard pattern
 Refractive index (n 20/D) : 1.400 - 1.404
 Density (at 25°C) : 0.868 - 0.878 g/mL
 Assay (GC) : min. 98.00 %

GHS Safety Information

Hazard Statement(s) : H226
 Signal Word : Warning
 Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 10-66
 S-Phrase(s) : 23-25
 WGK : 1
 RTECS : NS9800000
 Flash Point(°F) : 77°F
 Flash Point(°C) : 25°C
 Storage Temperature(°C) : Store below 30°C

Transport Information

UN No. : 1104
 Class : 3
 Packaging Group : 3
 Marine Pollutant : No



Product Information

ADR/RID : 1104 3/PG 3
IMDG : 1104 3/PG 3
IATA : 1104 3/PG 3

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein