
Eosin yellow, Practical grade (Spirit soluble, Free acid)

GRM1060

Product Identifier

CAS No.	:	15086-94-9
EC No.	:	239-138-3
C.I. No.	:	45380:2
Molecular Formula	:	C ₂₀ H ₈ Br ₄ O ₅
Molecular Weight	:	647.89
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Pink to orange to brown powder
Solubility	:	10 ppm solution in methanol yields clear orange yellow solution
Loss on drying (at 110°C, 1 hr)	:	<= 2.40%
Absorption maxima (10 ppm in methanol)	:	520 - 530 nm

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LM5800000
WGK	:	1