
Sodium malonate dibasic monohydrate, Bacteriological

GRM10683

Product Identifier

CAS No.	:	26522-85-0
EC No.	:	205-514-0
Molecular Formula	:	$C_3H_2Na_2O_4 \cdot H_2O$
Molecular Weight	:	166.04
HS Code	:	2917 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
Assay (NT)	:	99.00 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3