
Choline chloride, 98%

GRM10699


Product Identifier

CAS No.	:	67-48-1
EC No.	:	200-655-4
Molecular Formula	:	C ₅ H ₁₄ ClNO
Molecular Weight	:	139.62
HS Code	:	2923 10 00
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	White deliquescent crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 5.00%
Assay (NT/AT, on anhydrous basis)	:	98.00 - 101.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	KH2975000
WGK	:	1