

---

## Paraffin wax pellets

**GRM10702W**

---

### Product Identifier

CAS No.	:	8002-74-2
EC No.	:	232-315-6
HS Code	:	2712 20 00
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Colourless or white translucent pellets
Solubility	:	33.3 mg soluble in 1 mL of benzene
Congealing Point	:	58 - 60°C

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RV0350000
WGK	:	nwg