

RM1074 Phenylmercuric acetate

Product Number Packing

RM1074	: 100G
RM1074	: 25G

Product Information

Product Code	: RM1074
Product Name	: Phenylmercuric acetate
Synonym	: Mercury phenyl acetate; Phenylmercury acetate
Molecular Formula	: $C_8 H_8 H g O_2$
Molecular Weight	: 336.74
CAS No.	: 62-38-4
EC No.	: 200-532-5
HS Code	: 2852 10 00
Shelf Life	: 4 years

Technical Specification

Appearance	: White to yellow crystals or powder
Solubility	: 33.3 mg soluble in 1 mL of alcohol
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 148 - 151°C
Assay (KSCN Titration)	: min. 98.00 %

GHS Safety Information

Hazard Statement(s)	: H301-H314-H372-H410
Precautionary Statement(s)	: P273-P280-P301+P310-P305+P351+P338-P310-P501
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 25-34-48/24/25-50/53
S-Phrase(s)	: 23-24/25-37-45-60-61
WGK	: 3
RTECS	: OV6475000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	



Transport Information

UN No.	: 1674
Class	: 6.1



Product Information

Packaging Group	:	2
Marine Pollutant	:	YES
ADR/RID	:	1674 6.1/PG 2
IMDG	:	1674 6.1/PG 2
IATA	:	1674 6.1/PG 2

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein