

RM1082 Tungstic acid sodium salt dihydrate, Hi-AR™/ACS

Product Number Packing

RM1082	: 100G
RM1082	: 500G

Product Information

Product Code	: RM1082
Product Name	: Tungstic acid sodium salt dihydrate, Hi-AR™/ACS
Synonym	: Sodium tungstate dihydrate, Hi-AR™/ACS
Molecular Formula	: $N a_2 W O_4 \cdot 2 H_2 O$
Molecular Weight	: 329.85
CAS No.	: 10213-10-2
EC No.	: 236-743-4
HS Code	: 2841 80 00
Shelf Life	: 4 years

Technical Specification

Appearance	: Colourless or white to slightly yellow crystals or powder
Solubility	: 100 mg soluble in 1 mL of water
pH (2% in water at 25°C)	: 7.00 - 10.00
Chloride (Cl)	: $\leq 0.005\%$
Arsenic (As)	: $\leq 0.0005\%$
Lead (Pb)	: $\leq 0.001\%$
Iron (Fe)	: $\leq 0.001\%$
Copper (Cu)	: $\leq 0.001\%$
Nitrogen compounds (as N)	: $\leq 0.001\%$
Sulphate (SO ₄)	: $\leq 0.01\%$
Water insoluble matter	: $\leq 0.01\%$
Assay (Lead Nitrate Titration)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H302
Signal Word	: Warning
Hazard Pictogram(s)	: 

Risk and Safety Information

R-Phrase(s)	: 22
S-Phrase(s)	: 36
WGK	: 1
RTECS	: YO7900000
Storage Temperature(°C)	: Store below 30°C

Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods