

GRM1085**Sodium salicylate, Hi-AR™****Product Number****Packing**

GRM1085 : 250G

Product Information

Product Code : GRM1085
Product Name : Sodium salicylate, Hi-AR™
Synonym : 2-Hydroxybenzoic acid sodium salt; salicylic acid sodium salt
Molecular Formula : C₇H₅NaO₃
Molecular Weight : 160.10
CAS No. : 54-21-7
EC No. : 200-198-0
HS Code : 2918 21 00
Shelf Life : 4 years

Technical Specification

Appearance : Colourless or white to off-white crystals or powder or flakes
Solubility : 33.3 mg soluble in 1 mL of water
Chloride (Cl) : ≤ 0.02%
Sulfate (SO₄) : ≤ 0.06%
Loss on drying (at 105°C, 2 hr) : ≤ 0.50%
Assay (NT) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302-H319
Precautionary Statement(s) : P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-36
S-Phrase(s) : 26
WGK : 1
RTECS : VO5075000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods



Product Information

IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein