

Immersion oil, High viscosity

GRM10911

Product Identifier

HS Code : 3403 99 00
Storage : Below 30°C
Shelf life : 2 years

Technical Specification

Appearance : Colourless to pale greenish yellow viscous liquid
Solubility : 33.3 mg soluble in 1 mL of xylene
Refractive index (n 20/D) : 1.5110 - 1.5210
Viscosity : 1125 - 1375 cSt (Spindle no. : 18, Speed : 1.5 rpm, at 23°C)

Safety Information

Hazard Pictogram(s) :



Signal Word : Warning
Hazard Statement(s) : H315- H410
Precautionary Statement(s) : P273- P280- P302 + P352- P332 + P313- P362 + P364- P501
UN No. : 3082
Class : 9
Packing Group : III
RTECS : Not available
WGK : 2