

Potassium bromate

GRM1092


Product Identifier

CAS No.	:	7758-01-2
EC No.	:	231-829-8
Molecular Formula	:	KBrO ₃
Molecular Weight	:	167.00
HS Code	:	2829 90 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Iodimetry)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H271- H301- H350
Precautionary Statement(s)	:	P201- P210- P280- P301 + 310- P308 + P313- P370 + P378
UN No.	:	1484
Class	:	5.1
Packing Group	:	II
RTECS	:	EF8725000
WGK	:	3