
Bromoform, For mineral separation

GRM10947


Product Identifier

CAS No.	:	75-25-2
EC No.	:	200-854-6
Molecular Formula	:	CHBr ₃
Molecular Weight	:	252.73
HS Code	:	2903 39 19
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow liquid.
Solubility	:	1 mL miscible in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n _{20/D})	:	1.5900 - 1.6000
Density (at 25°C)	:	2.885 - 2.895 g/mL
Stabilizer	:	~1.00% ethanol
Assay (GC)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H315- H319- H331- H411
Precautionary Statement(s)	:	P261- P273- P305+P351+P338- P311
UN No.	:	2515
Class	:	6.1
Packing Group	:	III
RTECS	:	PB5600000