

GRM1100 Ammonium bromide, Hi-LR™**Product Number Packing**

GRM1100 : 500G

Product Information

Product Code : GRM1100
Product Name : Ammonium bromide, Hi-LR™
Molecular Formula : NH₄Br
Molecular Weight : 97.94
CAS No. : 12124-97-9
EC No. : 235-183-8
HS Code : 2827 59 00
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 4.80 - 6.00
Chloride (Cl) : ≤ 0.5%
Iron (Fe) : ≤ 0.001%
Sulfate (SO₄) : ≤ 0.05%
Assay (AT) : 99.00 - 101.00 %

GHS Safety Information

Hazard Statement(s) : H319
Precautionary Statement(s) : P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

RTECS : BO9155000
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods