
Micromountain

GRM11015

Product Identifier

HS Code : 3821 00 00
Storage : Below 30°C
Shelf life : 2 years

Technical Specification

Appearance : Colourless to pale greenish yellow viscous liquid
Solubility : 33.3 mg soluble in 1 mL of xylene
Density (at 23°C) : ~ 0.923 g/mL
Refractive index (n 20/D) : 1.5110 - 1.5210
Viscosity : max. 150 cSt

Safety Information

Hazard Pictogram(s) :



Signal Word : Warning
Hazard Statement(s) : H315- H319- H335
Precautionary Statement(s) : P261- P305+P351+P338
UN No. : 3082
Class : 9
Packing Group : III
RTECS : Not available
WGK : 3