
Benzoin

GRM1106

Product Identifier

CAS No.	:	119-53-9
EC No.	:	204-331-3
Molecular Formula	:	C ₁₄ H ₁₂ O ₂
Molecular Weight	:	212.25
Synonym	:	α -Hydroxy- α -phenylacetophenone
HS Code	:	2914 40 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light yellow crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of acetone
FTIR	:	Matches with the standard pattern
Melting range	:	129 - 138°C
Assay (GC)	:	min. 98.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DI1590000
WGK	:	2