

Sodium tetraphenylborate, Hi-AR™/ACS

GRM1156


Product Identifier

CAS No.	:	143-66-8
EC No.	:	205-605-5
Molecular Formula	:	C ₂₄ H ₂₀ BNa
Molecular Weight	:	342.22
Synonym	:	Kalignost, Hi-AR™/ACS
HS Code	:	2931 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to beige crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of acetic acid
Loss on drying (at 105°C, constant weight)	:	<= 0.5%
Clarity of solution	:	Passes test
Assay (NT/Gravimetry)	:	99.50 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301
Precautionary Statement(s)	:	P264- P270- P301 + P310- P330- P405- P501
UN No.	:	2811
Class	:	6.1
Packing Group	:	III
RTECS	:	ED3362500
WGK	:	1