
SPADNS trisodium salt, Hi-ARTM/ACS

GRM1157


Product Identifier

CAS No.	:	23647-14-5
EC No.	:	245-803-9
Molecular Formula	:	C ₁₆ H ₉ N ₂ Na ₃ O ₁₁ S ₃
Molecular Weight	:	570.41
Synonym	:	Sulphanilic acid azochromotrop, Hi-AR TM /ACS
HS Code	:	2927 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Red to brown powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Clarity	:	Passes test
Sensitivity to fluoride	:	Passes test

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1