

## Tetrabutylammonium hydrogensulfate, Hi-LR<sup>®</sup>

**GRM1163**


### Product Identifier

CAS No.	:	32503-27-8
EC No.	:	251-068-5
Molecular Formula	:	C <sub>16</sub> H <sub>37</sub> NO <sub>4</sub> S
Molecular Weight	:	339.53
Synonym	:	Tetrabutylammonium bisulfate
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless to off-white crystals or powder or lumps
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration)	:	97.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H314- H412
Precautionary Statement(s)	:	P260- P273- P280- P301+P330+P331- P305+P351+P338- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	2