
Tetrabutylammonium iodide, Hi-AR™

GRM1164


Product Identifier

CAS No.	:	311-28-4
EC No.	:	206-220-5
Molecular Formula	:	C ₁₆ H ₃₆ IN
Molecular Weight	:	369.37
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to tan crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of hot water
FTIR	:	Matches with the standard pattern
Melting range	:	141 - 150°C
Assay (AT)	:	99.00 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319
Precautionary Statement(s)	:	P270- P280- P301 + P312- P302 + P352- P305+P351+P338- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	BS5450000
WGK	:	2