

---

## Sodium pyruvate, Hi-AR™

**GRM1181WK**

---

### Product Identifier

CAS No.	:	113-24-6
EC No.	:	204-024-4
Molecular Formula	:	C <sub>3</sub> H <sub>3</sub> NaO <sub>3</sub>
Molecular Weight	:	110.04
Synonym	:	Pyruvic acid sodium salt, Hi-AR™
HS Code	:	2918 30 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	Yellowish white, homogeneous, free flowing powder.
Solubility	:	33.3 mg soluble in 1 mL of water
Clarity	:	1% w/v aqueous solution is clear without any haziness.
pH (10% in water at 25°C)	:	5.50 - 7.00
FTIR	:	Matches with the standard pattern
Loss on drying (at 105°C, 2 hr)	:	<= 0.5%
Assay (NT, on dry basis)	:	min. 98.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1