

## Sodium hydroxide pellets, Hi-LR<sup>TM</sup>

**GRM1183**

### Product Identifier

CAS No.	:	1310-73-2
EC No.	:	215-185-5
EC Index No	:	011-002-00-6
Molecular Formula	:	NaOH
Molecular Weight	:	40.0
HS Code	:	2815 11 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White hygroscopic pellets
Solubility	:	100 mg soluble in 1 mL of water
Carbonate (Na <sub>2</sub> CO <sub>3</sub> )	:	<= 2.0%
Chloride (Cl)	:	<= 0.02%
Iron (Fe)	:	<= 0.005%
Sulfate (SO <sub>4</sub> )	:	<= 0.01%
Assay (HCl Titration)	:	97.00 - 103.50%

### Safety Information

Hazard Pictogram(s) :



Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P305+P351+P338- P310
UN No.	:	1823
Class	:	8
Packing Group	:	II
RTECS	:	WB4900000
WGK	:	1